

Deep Lab Cut

INSTALL

[Install Deep Lab Cut 2.0](#) :

Preferably on an Ubuntu system (the plan is to train the model, get the weights, and then transfer the weights to a Docker which is based on an Ubuntu OS. It's easier to transfer the weights if both are done on Ubuntu systems)

COMPLETE GUIDE

Refer: <https://github.com/AlexEMG/DeepLabCut/blob/master/docs/functionDetails.md>

Example notebooks: <https://github.com/AlexEMG/DeepLabCut/tree/master/examples>

TRAIN

- Train the model following instructions in the notebook `Deeplabcut_new.ipynb` in the git repo (https://github.com/ambareeshsrja16/Surgical_Tool_Tracking/blob/master/Deeplabcut_new.ipynb)
- Project config file is stored in the `path_config_file` variable - after you make the project, print this variable and store it, so that you can reuse it again if need be.
- For changing config.yaml file according to specific needs refer to <https://github.com/AlexEMG/DeepLabCut/blob/master/docs/functionDetails.md> (Read BOX 1)

```
# Training,Evaluation and Analysis configuration
```

```
TrainingFraction :- 0.95    # how much of the labelled data needs to be used for training
```

```
iteration: 5    #New weights will be saved with this number in the name
```

- Before you run this command - `deeplabcut.train_network(path_config_file)` - make sure to change some parameters in `pose_cfg.yaml` at the following path:

```
C:\....\DaVinci-Ambar-2020-02-14\dlc-models\iteration-2\DaVinciFeb14-trainset40shuffle1\train
```

Change these two: `all_joints`, `all_joints_names`; according to the robot.

```
all_joints:
```

```
-- 0
```

```
-- 1
-- 2
-- 3
-- 4
-- 5
-- 6
-- 7
-- 8
-- 9
-- 10
-- 11
all_joints_names:
- PSM1-roll_2
- PSM1-pitch_1
- PSM1-pitch_2
- PSM1-pitch_3
- PSM1-yaw_1
- PSM1-yaw_2
- PSM2-roll_2
- PSM2-pitch_1
- PSM2-pitch_2
- PSM2-pitch_3
- PSM2-yaw_1
- PSM2-yaw_2
```

(The numbers show the order in which you'll label these parts later)

Change these two parameters to monitor the training:

display_iters: 100 (so that it prints loss values after every 100 iterations, and you can track the progress)

save_iters: 500 (so that the model can be saved after every 500 iterations or so)

- For CUDA errors that pop up, use the following:

```
import tensorflow as tf
gpu_options = tf.GPUOptions(per_process_gpu_memory_fraction=0.6)
sess = tf.Session(config=tf.ConfigProto(gpu_options=gpu_options))
```

The above prevents TensorFlow from using up the entire GPU Memory

EVALUATE

Models can be evaluated using DLC's inbuilt evaluation functions. Some additional graphs created for the SuperDeep paper were made using functions in the `Deeplabcut_new.ipynb` notebook

DOCKER

After training the model, you should have some file structure like this:

PC > OS (C:) > Users > asree > PycharmProjects > DEEPLABCUT_New > DaVinci-Ambar-2020-02-14				
Name	Date modified	Type	Size	
.ipynb_checkpoints	2/26/2020 12:59 PM	File folder		
data_Menglong_Ye	2/21/2020 10:43 AM	File folder		
dlc-models	2/27/2020 11:12 PM	File folder		
labeled-data	2/28/2020 7:07 PM	File folder		
training-datasets	2/27/2020 11:12 PM	File folder		
videos	2/14/2020 4:18 PM	File folder		
backup_txt	2/27/2020 5:35 PM	TXT File	7 KB	
config.yaml	2/28/2020 7:26 PM	YAML File	2 KB	
Joint_Angles_List.pkl	2/21/2020 2:47 PM	PKL File	45 KB	
Pose_to_Joint_Angles_notebook	2/27/2020 2:10 PM	IPYNB File	140 KB	
PSM1_joint_angles.npy	2/26/2020 1:15 PM	NPY File	21 KB	
PSM2_joint_angles.npy	2/26/2020 1:16 PM	NPY File	21 KB	
sample	8/22/2017 2:48 AM	XML Document	483 KB	
sim-PSM1-2016Feb05T175836	8/22/2017 2:48 AM	XML Document	483 KB	
sim-PSM2-2016Feb05T175836	8/22/2017 2:48 AM	XML Document	482 KB	
Utility	2/27/2020 2:10 PM	IPYNB File	160 KB	

The project folder has to be copied to Docker for running the prediction, after the Docker container is up and running!

(The folder `dlc-models` has all the weights, numbered according to the iteration (you can change this in the `config.yaml` file when you are training). But we copy the entire project folder to make the other steps below easier. You can obviously change this when you understand how the whole thing works)

Clone the git repo : https://github.com/ambareeshsija16/Surgical_Tool_Tracking

Go to (cd to) the folder containing Dockerfile, and run the following commands

1. `docker image build -t dlc_username/dlcdocker .`
builds the docker container from image
2. `docker run -it --gpus "device=2" dlc_username/dlcdocker bash`
runs the docker container
#If you want to use all GPUs use:
`docker run -it --gpus 'all' dlc_username/dlcdocker bash`

Copy the project folder to Docker using:

```
docker cp test-jingpei-2020-01-09  
elegant_albattani:root/DLCROS_ws/src/Surgical_Tool_Tracking/ForwardPassDeepLabCut
```

Note:

test-jingpei-2020-01-09 - name of the project folder (in the picture above the name is *DaVinci-Ambar-2020-02-14*)
elegant_albattani - name of your docker container

Run these commands (they are hacks to prevent some errors - either in DLC, or the ROS2-Python3 binding - hopefully in the future these will clear up. For more details refer the last parts of this document, those are not formatted though)

```
export DLClight=True  
rm -rf /opt/ros/kinetic/lib/python2.7/dist-packages/cv2.so
```

For running the evaluation/ prediction/ forward pass code

```
cd ~/DLCROS_ws/Surgical_Tool_Tracking/ForwardPassDeepLabCut  
python3 generator_clipped_image...
```

COMMON ERRORS

- “FileNotFoundError: It seems the model for shuffle s and train ...”
Check the config.yaml file that you copied inside the Docker, and check the **project_path** in this .yaml file, it could be incorrect
- “No module name numpy.testing.decorators”
Uninstall numpy version 1.18.0 and install

```
pip uninstall numpy  
pip install numpy==1.17.0
```


<https://blog.csdn.net/AiBigData/article/details/103738505>

DLC : MISCELLANEOUS

The following is a rough clipboard with information that was used during the project

Test Docker

```
sudo docker run hello-world
```

Test NVIDIA, Jacob's command

```
sudo docker run --gpus all nvidia/cuda nvidia-smi
```

Added my username to the docker group:

```
sudo usermod -aG docker ambareesh
```

Usually, you have to restart, to make sure the changes are reflected, but check on a new terminal if you are able to run the hello-world docker test without sudo, if yes, you are good

Final

Ran 1, 2

1)

```
docker run --gpus all -t -d -p 2351:8888 -e USER_HOME=$HOME/Docker4DeepLabCut2.0  
dlc_username/dlcdocker bash
```

-d = detach, run container in bg, and print id

-p = attach the port

-e = set environment variables

-t = allocate a pseudo tty (Exits without this for some reason)

dlc_username/dlcdocker = image

bash -> add at the end to keep the container running!

Serverhost:2351 linked to container:8888

2)

```
docker exec -it sweet_burnell /bin/bash
```

-i= interactive

clever_rosalind= random name, (if you use --user name when you start the container you can choose your own)

`docker exec --user ambareesh -it sweet_burnell /bin/bash`
unable to find user ambareesh: no matching entries in passwd file?

3) `cd ../../../../home/`
4) `export DLClight=True`
5) `ipython -> import deeplabcut`

6) For Jupyter notebook

Connect your localhost:2351 to the serverhost:2351

`ssh -L localhost:2351:localhost:2351 ambareesh@unicorn2.ucsd.edu`

Run `jupyter-notebook --allow-root` from inside container

You'll get a token (but with 8888, this is what the container is generating, this is routed to 2351 of server host and in turn to your 2351, so copy-paste the token to your browser and replace 8888 with 2351)

Other

`docker ps` (show only running ones) `-a` (show everything)
`docker container prune`
`docker stop <container_name>`

Accessing the Server

`ambareesh@unicorn2.ucsd.edu`
perception

Image building

```
docker image build -t dlc_username/dlcdocker .
```

Running the docker

- `docker run -it dlc_username/dlcdocker`
- `docker run -it --gpus "device=2" dlc_username/dlcdocker bash`
- `docker run -it --gpus 'all' dlc_username/dlcdocker bash`

Do first if you have added CMD["bash"] in the docker file!, if not then do the second one

Copy stuff from docker to desktop

Docker cp <name>:~/path host_path

ROS

To show what is being published

Rosrun image_view image_view image:=/dlc_prediction_topic

To record a view from the bag

Rosrun image_view video_recorder _filename:="/home/.._video.avi" image:=(topic)

Copy files from ssh to my system

scp ambareesh@unicorn2.ucsd.edu:~/version_1_Docker/Docker4DeepLabCut2.0/Dockerfile .

scp ambareesh@unicorn2.ucsd.edu:~/version_1_Docker/Docker4DeepLabCut2.0/Dockerfile .
scp

Users/ambareeshsnjayakumari/Desktop/Surgical_Robot/Git/Surgical_Tool_Tracking/ForwardPassDeepLabCut/ ambareesh@unicorn2.ucsd.edu:~/version_1_Docker/Docker4DeepLabCut2.0/

Method 1

1. Remove last two lines for sourcing the ros stuff from the docker file, and then build it
Now run using gpus, export DLClight=True
Try python3, import deeplabcut (should work with np warnings), and then import cv_bridge
(wouldn't work),
Import sys
sys.path.append("/opt/ros/kinetic/lib/python2.7/dist-packages") now import cv_bridge and it
should work. When opening new terminals to make the ROS stuff work make sure that you
source the stuff correctly.

2. Add the last two lines in the dockerfile as such and then remove cv2.so from the
/opt/ros/kinetic/lib/python2.7/dist-packages folder, this might work.

TO RUN

scp weights to host

docker cp

/home/jingpei/test-jingpei-2020-01-09/dlc-models/iteration-0/testJan9-trainset95shuffle1/train/snapshot-7100.meta elegant_albattani:root/weights/snapshot-7100.meta

```
docker image build -t dlc_username/dlcdocker .  
docker run -it --gpus "device=2" dlc_username/dlcdocker bash
```

```
export DLClight=True  
rm -rf /opt/ros/kinetic/lib/python2.7/dist-packages/cv2.so  
cd ~/DLCROS_ws/Surgical_Tool_Tracking/ForwardPassDeepLabCut
```

```
python3 prediction_v1.py
```

```
cp prediction_v1.py prediction_v2.py  
vim prediction_v2.py
```

```
"/root/DLCROS_ws/Surgical_Tool_Tracking/ForwardPassDeepLabCut/DaVinci-Ambar-2019-10-31/config.yaml"
```

```
vim /usr/local/lib/python3.6/dist-packages/deeplabcut/utils/auxiliaryfunctions.py  
Use four spaces
```

```
python3 prediction_v2.py
```

```
vim  
~/DLCROS_ws/Surgical_Tool_Tracking/ForwardPassDeepLabCut/DaVinci-Ambar-2019-10-31/  
config.yaml
```

```
project_path  
/root/DLCROS_ws/Surgical_Tool_Tracking/ForwardPassDeepLabCut/DaVinci-Ambar-2019-10-31
```

```
video_sets  
/root/DLCROS_ws/Surgical_Tool_Tracking/ForwardPassDeepLabCut/DaVinci-Ambar-2019-10-31/videos/video_1.avi
```

Snapshot index problem

```
~/DLCROS_ws/Surgical_Tool_Tracking/ForwardPassDeepLabCut/DaVinci-Ambar-2019-10-31/  
dlc-models/iteration-2/DaVinciOct31-trainset95shuffle1/train  
Added this to dockerfile, scp ed the file to host, and copied it in.
```


io.imread problem for PNG files, try the below thing
image = io.imread(image_name, plugin='matplotlib')

ssh ambareesh@unicorn2.ucsd.edu

DLC

Report error in

ruamel.yaml.scanner.ScannerError: mapping values are not allowed here
in

"/root/DLCROS_ws/Surgical_Tool_Tracking/ForwardPassDeepLabCut/DaVinci-Ambar-2019-10-31/config.yaml", line 9, column 116

"/usr/local/lib/python3.6/dist-packages/deeplabcut/utils/auxiliaryfunctions.py", line 120, in
read_config

cfg = ruamelFile.load(f)

Need to write an else statement for other errors

MULTIPLE TERMINALS

You can run `docker exec -it <container> bash` from multiple terminals to launch several sessions connected to the same container.

No module name numpy.testing.decorators

Uninstall numpy version 1.18.0 and install

pip uninstall numpy

pip install numpy == 1.17.0

<https://blog.csdn.net/AiBigData/article/details/103738505>

Running on a fresh system

- Download/clone from git https://github.com/ambareeshsrja16/Surgical_Tool_Tracking
- Run
 - docker image build -t dlc_username/dlcdocker .
 - docker run -it dlc_username/dlcdocker

OR

```
docker run -it --gpus "device=2" dlc_username/dlcdocker bash
OR
```

```
docker run -it --gpus 'all' dlc_username/dlcdocker bash
```

- Will face this error: *Python 3.6 PPA removed from public access*
Solve using: <https://github.com/DeepLabCut/Docker4DeepLabCut2.0/issues/21>
- If intel-openmp distribution error occurs do this:
<https://github.com/DeepLabCut/Docker4DeepLabCut2.0/issues/25>
- Copy test folder (test-Jingpe....) to the docker inside ForwardPassDeepLabcut
- Try Python3 generator_clipped_image...
 - Import cv2 error: remove cv2.so from the
/opt/ros/kinetic/lib/python2.7/dist-packages folder
 - No module named "wx"
export DLClight=True
 - No module name numpy.testing.decorators
Uninstall numpy version 1.18.0 and install
pip uninstall numpy
pip install numpy==1.17.0
<https://blog.csdn.net/AiBigData/article/details/103738505>
 - export DLClight=True
 - Copy test JingPei folder

Copy Weights to Docker!

```
docker cp test-jingpei-2020-01-09
```

```
elegant_albattani:root/DLCROS_ws/src/Surgical_Tool_Tracking/ForwardPassDeepLabCut
```